

**CONCEPTION TECHNOLOGIES
ENHANCE WATER FOR TISSUE CULTURE**

Catalog # WTC-1500 (1.5L)
FOR LABORATORY PROCEDURES ONLY
Other uses must be qualified by the end user

INTENDED USE

Water for Tissue Culture is intended for use in laboratory procedures requiring a non-pyrogenic, high purity grade of water. Water for Tissue Culture is produced from USP sterile water, which meets "Water for Injection" criteria. It has been found to exhibit no cytotoxic effects when tested by the mouse zygote QC test of Quinn et. al. (Fertil Steril 41:202, 1984).

COMPONENTS

Distilled Water, sterile filtered.
See Certificate of Analysis for specifications of individual lots.
Contains no bacteriostat, antimicrobial agent or added buffer.

STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers at room temperature, (below 30°C) for up to one year from date of manufacture. The expiration date is shown on the product label and Certificate of Analysis.

DIRECTIONS FOR USE

Each laboratory should make its own determination of suitability for use in each particular procedure.

NOT INTENDED FOR INJECTION BY USUAL PARENTERAL ROUTES IN HUMANS OR ANIMALS

PRECAUTIONS AND WARNINGS

Do not use any bottle of water showing evidence of particulate matter or cloudiness. To avoid problems with contamination, practice aseptic technique and discard minimal amounts of excess water remaining in the bottle.

RELATED PRODUCTS

Conception Technologies has a full line of products for the Reproductive Medical Specialist. Please call or write for specific information or to receive a copy of our current price list. For technical questions, or to reach our Customer Service Department:



**CONCEPTION
TECHNOLOGIES**

858.824.0888/888.995.8081
6835 Flanders Drive, Ste 500
San Diego, CA 92121
www.conceptiontechnologies.com